



Anti-KIF5C polyclonal antibody (CABT-BL2062)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human KIF5C.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Kinesin is a microtubule-associated force-producing protein that may play a role in organelle transport.
---------------------	--

GENE INFORMATION

Entrez Gene ID [3800](#)

Protein Refseq [NP_004513](#)

UniProt ID [O60282](#)
