



Anti-KIF5C polyclonal antibody (CABT-BL2063)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant GST fusion protein corresponding to 82 mouse Kif5c.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB, IHC
Format	Liquid
Buffer	In PBS (50% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Kinesin is a microtubule-associated force-producing protein that may play a role in organelle transport.
--------------	--

GENE INFORMATION

Entrez Gene ID	3800
Protein Refseq	NP_004513

UniProt ID

[O60282](#)
