



# Anti-KLHDC9 (aa 215-264) polyclonal antibody (CABT-BL2071)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide corresponding to a region within internal sequence amino acids 215-264 (AEPEVAGHWS HGKIKEPPV APHLMEQLAR LVSSGQGSQK GPHGLRHHSC) of Human KLHDC9.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Format</b>	Liquid
<b>Concentration</b>	Lot specific
<b>Buffer</b>	2% Sucrose, PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

## GENE INFORMATION

Entrez Gene ID [126823](#)

Protein Refseq [NP\\_001007256](#)

---

UniProt ID [Q8NEP7](#)

---