



Anti-KLHL24 (aa 190-220) polyclonal antibody (CABT-BL2073)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 190-220 of Human KLHL24, conjugated to KLH.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Perikaryon. Cell projection, axon. Cytoplasm
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Entrez Gene ID	54800
Protein Refseq	NP_060114

UniProt ID

[Q6TFL4](#)
