



Anti-KLHL28 (aa 362-411) polyclonal antibody (CABT-BL2074)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A region within synthetic peptide: DTNTWTSLER MNESRSTLGV VVLAGELYAL GGYDQGQSYLQ SVEKYIPKIR, corresponding to internal sequence amino acids 362-411 of Human BTB domain containing 5
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB, ELISA
Cellular Localization	Nuclear
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	KLHL28 (kelch-like family member 28) is a protein-coding gene. An important paralog of this gene is KLHL17.
---------------------	---

GENE INFORMATION

Entrez Gene ID [54813](#)

Protein Refseq [NP_060128](#)

UniProt ID [Q9NXS3](#)
