



Anti-KLHL31 (aa 501-550) polyclonal antibody (CABT-BL2076)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A region within synthetic peptide: TPRGWHCAVT LSDRVYVMGG SQLGPRGERV DVLTVCEYSP ATGQWSYAAP, corresponding to amino acids 501-550 of Human KLHL31
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	401265
Protein Refseq	NP_001003760
UniProt ID	Q9H511