



Anti-KLHL35 (aa 72-121) polyclonal antibody (CABT-BL2078)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 72-121(PRRFMDLAEV IVVIGGCDRK GLLKLPFADA YHPESQRWTP LPSLPGYTRS) of Human KLHL35.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	283212
Protein Refseq	NP_001034637

UniProt ID

[Q6PF15](#)
