



Anti-KRCC1 (aa 35-84) polyclonal antibody (CABT-BL2098)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N-terminal amino acids 35-84 (KMKGDYLETCGYKGEVNSRPTYRMFDQRLPSETIQTYPR SCNIPQTVENR) of Human KRCC1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	51315
Protein Refseq	NP_057702

UniProt ID

[A0A024R5P4](#)
