



Anti-Soybean KTi3 polyclonal antibody (CABT-BL2103)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Full length protein
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Soy Bean
Purification	IgG fraction
Conjugate	Unconjugated
Applications	WB, Dot, ELISA, Conjugation
Format	Liquid
Size	100 µg
Buffer	0.15M Sodium chloride, 0.02M Potassium phosphate. pH 7.2
Preservative	0.01% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Trypsin inhibitors are chemicals that reduce the availability of biologically active trypsin. Soybeans contain several inhibitors; the one in the chart is considered the primary one. All of them bind chymotrypsin to a lesser degree.
---------------------	--

GENE INFORMATION

Entrez Gene ID [547831](#)

Protein Refseq [NP_001238611](#)

UniProt ID [P01070](#)
