



Anti-L3MBTL2 (aa 613-662) polyclonal antibody (CABT-BL2108)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: LKAKEATKKK KKQFGKKRKR IPPTKTRPLR QGSKKPLLED DPQGARKISS, corresponding to amino acids 613-662 of Human L3MBTL2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Cellular Localization	Nucleus.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Lethal(3)malignant brain tumor-like 2 protein is a protein that in humans is encoded by the L3MBTL2 gene.
---------------------	---

GENE INFORMATION

Entrez Gene ID [83746](#)

Protein Refseq [NP_113676](#)

UniProt ID [Q969R5](#)
