



Anti-LCA5L (aa 621-700) polyclonal antibody (CABT-BL2118)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal amino acids 621 - 700 (GSEEPLQSKE SHPLPPSQAS TSHAFGDSKV TVVNSIKPSS PTEGKRKIII) of Human C21ORF13
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	150082
Protein Refseq	NP_689718

UniProt ID

[O95447](#)
