



Anti-LCN10 (aa 101-131) polyclonal antibody (CABT-BL2122)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within internal sequence amino acids 101-131 of Human LCN10
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Secreted.
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

BACKGROUND

Introduction	Members of the lipocalin family, such as LCN10, have a common structure consisting of an 8-stranded antiparallel;beta-barrel that forms a cup-shaped ligand-binding pocket or calyx. Lipocalins generally bind small hydrophobic;ligands and transport them to specific cells (Suzuki
---------------------	---

GENE INFORMATION

Entrez Gene ID	414332
----------------	------------------------

Protein Refseq	NP_001001712
----------------	------------------------------

UniProt ID	Q6JVE6
------------	------------------------
