



Anti-LCN6 (aa 108-157) polyclonal antibody (CABT-BL2133)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within the internal sequence amino acids 108-157 (WVLATNFRDY AIIFTQLEFG DEPFNTVELY SLTETASQEA MGLFTKWSRS) of human LCN6.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Secreted
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	158062
Protein Refseq	NP_945184

UniProt ID

[A0A024R8I9](#)
