



Anti-LHPP polyclonal antibody (CABT-BL2145)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|---------------------------|---|
| Immunogen | Recombinant protein corresponding to amino acids of human LHPP. |
| Isotype | IgG |
| Source/Host | Rabbit |
| Species Reactivity | Human, Rat |
| Purification | Antigen affinity purification |
| Conjugate | Unconjugated |
| Applications | WB, IHC, IF |
| Format | Liquid |
| Buffer | In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide) |
| Preservative | 0.02% Sodium Azide |
| Storage | Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing. |

BACKGROUND

| | |
|---------------------|--|
| Introduction | LHPP is a member of the HAD-like hydrolase superfamily. It is a unique enzyme that hydrolyzes not only oxygen-phosphorus bonds in inorganic pyrophosphate but also nitrogen-phosphorus bonds in phospholysine, phosphohistidine and imidodiphosphate in vitro. |
|---------------------|--|

GENE INFORMATION

Entrez Gene ID [64077](#)

Protein Refseq [NP_001161352](#)

UniProt ID [Q9H008](#)
