



Anti-LMNB1 polyclonal antibody (CABT-BL2156)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to C-terminal of human Lamin B1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Cellular Localization	Nucleus inner membrane.
Format	Liquid
Buffer	1% BSA, 50mM Tris, pH 7.6
Preservative	0.1% Sodium Azide
Storage	Store at +4°C.

BACKGROUND

Introduction	Lamins are components of the nuclear lamina, a fibrous layer on the nucleoplasmic side of the inner nuclear membrane, which is thought to provide a framework for the nuclear envelope and may also interact with chromatin.
--------------	--

GENE INFORMATION

Entrez Gene ID [4001](#)

Protein Refseq [NP_001185486](#)

UniProt ID [B4DZT3](#)
