



Anti-LMNB1 (C-terminal) polyclonal antibody (CABT-BL2159)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to C-terminus of human lamin B.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	IHC-P, WB
Cellular Localization	Nucleus Envelope
Preservative	None
Storage	2-8°C

BACKGROUND

Introduction	Lamins are components of the nuclear lamina, a fibrous layer on the nucleoplasmic side of the inner nuclear membrane, which is thought to provide a framework for the nuclear envelope and may also interact with chromatin.
--------------	--

GENE INFORMATION

Entrez Gene ID	4001
----------------	----------------------

Protein Refseq [NP_001185486](#)

UniProt ID [B4DZT3](#)
