



Anti-LNP1 (internal region) polyclonal antibody (CABT-BL2165)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Peptide with sequence SSDCHPRRSHEDQ, from the internal region of human LNP1 according to NP_001078920.1.
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	0.5% BSA, Tris buffered saline, pH 7.3
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	LNP1 (leukemia NUP98 fusion partner 1) is a protein-coding gene. Diseases associated with LNP1 include leukemia, and acute lymphoblastic leukemia.
---------------------	--

GENE INFORMATION

Entrez Gene ID [348801](#)

Protein Refseq [NP_001078920](#)

UniProt ID [A1A4G5](#)
