



Anti-LOC391763 (full length) polyclonal antibody (CABT-BL2170)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Full length protein (Human, aa 1-198)
Isotype	IgG
Source/Host	Mouse
Species Reactivity	Human
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [391763](#)