



Anti-LOC396260 polyclonal antibody (CABT-BL2176)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Avidin
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Chicken
Purification	IgG fraction
Conjugate	Unconjugated
Applications	WB, IP, ELISA, Conjugation, Dot
Cellular Localization	Secreted
Format	Liquid
Buffer	0.15M Sodium chloride, 0.02M Potassium phosphate, pH 7.2
Preservative	0.01% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	396260
Protein Refseq	NP_990651

UniProt ID

[P02701](#)
