



Anti-LOC396260 (native) polyclonal antibody (CABT-BL2177)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Native purified Avidin from hen egg white.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Chicken
Conjugate	Unconjugated
Applications	WB, IHC, ELISA
Format	Liquid
Buffer	In 20 mM potassium phosphate buffer, 150 mM NaCl, pH 7.2 (0.01% sodium azide)
Preservative	0.01% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Entrez Gene ID	396260
Protein Refseq	NP_990651.1
UniProt ID	P02701