



Anti-LOC647589 (aa 51-100) polyclonal antibody (CABT-BL2180)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 51-100 (LRLHLLDNAD VALACARVLD QQEQQAACR LLEGCVPGG SQELVQLWND) of Human LOC647589
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	647589
Protein Refseq	NP_001177983

