



Anti-LOC652119 (aa 362-411) polyclonal antibody (CABT-BL2184)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within the C terminal amino acids 362-411 (RSSALPSGLL LHELLVSQEF LQQGLPFLET EAPGELEALE EAASLEEPLS) of Human LOC652119 (XP_946548).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [652119](#)