



Anti-LOC653186 (aa 1-50) polyclonal antibody (CABT-BL2185)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide: MKVINAKLDG FYSFPIFLFQ FLQATAQEEG IFECADPKLA ISAIWTFRDL, corresponding to amino acids 1-50 of Human LOC653186
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB, ELISA
Reconstitution	Add 100 µl of distilled water for a final concentration of 1mg/ml.
Format	Lyophilised
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [653186](#)