



Anti-LOC727787 (aa 271-333) polyclonal antibody (CABT-BL2186)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within amino acids 271 and 333 of Human LOC727787 (XP_001718677).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ICC/IF
Format	Liquid
Buffer	10% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7
Preservative	None
Storage	Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

GENE INFORMATION

Entrez Gene ID	727787
----------------	------------------------