



# Anti-LOC729020 (aa 187-215) polyclonal antibody (CABT-BL2188)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide conjugated to KLH, from within C terminal amino acids 187-215 of Human hCG_2024410 (Q2QD12)
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Mouse, Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Format</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Store at 4°C (up to 6 months). For long term storage store at -20°C

## GENE INFORMATION

<b>Entrez Gene ID</b>	<a href="#">729020</a>
<b>Protein Refseq</b>	<a href="#">NP_001137381</a>
<b>UniProt ID</b>	<a href="#">Q2QD12</a>