



Anti-LOQS polyclonal antibody (CABT-BL2192)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fragment derived from within residues 300 to the C-terminus of Fruit fly (<i>Drosophila melanogaster</i>) Loquacious.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Fruit fly
Purification	Protein A purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Nuclear
Format	Liquid
Buffer	1% BSA, PBS, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	34751
Protein Refseq	NP_001033903

UniProt ID

[Q4T2M6](#)
