



Anti-LRPAP1 polyclonal antibody (CABT-BL2196)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein corresponding to amino acids of human LRPAP1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC, IF
Format	Liquid
Buffer	In PBS, pH 7.5 (40% glycerol, 0.02% sodium azide)
Preservative	0.02% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	Interacts with LRP1/alpha-2-macroglobulin receptor and glycoprotein 330.
--------------	--

GENE INFORMATION

Entrez Gene ID [4043](#)

Protein Refseq [NP_002328](#)

UniProt ID [P30533](#)
