



Anti-LRRC25 (aa 36-85) polyclonal antibody (CABT-BL2208)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 36-85 (AEFSATCLNF SGLSLSLPHN QSLRASNVIL LDLSGNGLRE LPVTFFAHLQ) of Human LRRC25.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane; Single-pass type I membrane protein.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [126364](#)

Protein Refseq [NP_660299](#)

UniProt ID [Q8N386](#)
