



Anti-LRRC37A3 (internal region) polyclonal antibody (CABT-BL2218)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within the internal amino acids 1080-1129 (YTSTELIIEP EEPDSSGIN LSGFGSEQLD TNDESDVTST LSYILPYFSA) of Human LRRC37A3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane; single pass type 1 membrane protein
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [374819](#)

Protein Refseq [NP_955372](#)

UniProt ID [O60309](#)
