



Anti-LRRC4B (aa 88-118) polyclonal antibody (CABT-BL2219)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, between amino acids 88-118 from the N terminal region of Human LRRC4B.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane. Cell junction; synapse; presynaptic cell membrane.
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Entrez Gene ID	94030
Protein Refseq	NP_001073926

UniProt ID

[A0A024R4I8](#)
