



Anti-LRRC52 (N-terminal) polyclonal antibody (CABT-BL2220)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N-terminal amino acids QEVICTGKQLTEYPLDIPLNTRRLFLNENRITSLPAMHLG LLSDLVYLD Cof Human LRRC52.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	440699
Protein Refseq	NP_001005214

UniProt ID

[Q8N7C0](#)
