



Anti-LRRC63 (aa 491-521) polyclonal antibody (CABT-BL2222)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A KLH conjugated synthetic peptide corresponding to a region between C terminal amino acids 491-521 of human LRRC63. (Q05C16)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	220416
UniProt ID	J3KQU2