



Anti-LRRC66 (aa 756-805) polyclonal antibody (CABT-BL2223)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 756-805 (NVTFQTIPGK CKNQEDPFEK PLISAPDSGM YKTHLENASD TDRSEGLSPW) of Human LRRC66.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Membrane; Multi-pass membrane protein.
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID [339977](#)

Protein Refseq [NP_001019782](#)

UniProt ID [Q68CR7](#)
