



Anti-LRRC73 (aa 180-229) polyclonal antibody (CABT-BL2227)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within internal amino acids 180 - 229 (LDYNPLGDHV AGMLAVAVAS SRTLEVLDLE GTGLTNQSAQ TLLDMVENYP) of Human C6ORF154.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	221424
Protein Refseq	NP_001012992

UniProt ID

[Q5JTD7](#)
