



## Anti-ltp-BR1 (aa 1-15) polyclonal antibody (CABT-BL2243)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Immunogen	Synthetic peptide, corresponding to amino acids 1-15 of Barley (Hordeum vulgare) Lipid transfer protein
Isotype	IgG1
Source/Host	Goat
Species Reactivity	Barley
Purification	Protein A purified
Conjugate	Unconjugated
Applications	ELISA
Format	Liquid
Buffer	Tris buffered saline, pH 7.4
Preservative	None
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

### GENE INFORMATION

Entrez Gene ID [548228](#)