



Anti-LYPD6 (aa 130-160) polyclonal antibody (CABT-BL2245)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within C terminal amino acids 130-160 of Human LYPD6.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Cellular Localization	Secreted
Format	Liquid
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

BACKGROUND

Introduction	Members of the LY6 protein family (see SLURP1; MIM 606119), such as LYPD6, have at least one 80-amino acid LU domain; that contains 10 conserved cysteines with a defined disulfide-bonding pattern (Zhang et al., 2010 (PubMed:19653121)).
---------------------	---

GENE INFORMATION

Entrez Gene ID [130574](#)

Protein Refseq [NP_919298](#)

UniProt ID [Q86Y78](#)
