



# Anti-LYPD6 (aa 108-157) polyclonal antibody (CABT-BL2246)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide corresponding to a region within internal sequence amino acids 108 - 157 (RDSEHEGHKV CTSCCEGNIC NLPLPRNETD ATFATTSPIN QTNGHPRCMS) of human LYPD6.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Cellular Localization</b>	Secreted
<b>Format</b>	Liquid
<b>Buffer</b>	2% Sucrose, PBS
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

## BACKGROUND

<b>Introduction</b>	Members of the LY6 protein family (see SLURP1; MIM 606119), such as LYPD6, have at least one 80-amino acid LU domain;that contains 10 conserved cysteines with a defined disulfide-
---------------------	---

## GENE INFORMATION

Entrez Gene ID	<a href="#">130574</a>
----------------	------------------------

---

Protein Refseq	<a href="#">NP_919298</a>
----------------	---------------------------

---

UniProt ID	<a href="#">Q86Y78</a>
------------	------------------------

---