



# Anti-M (internal region) polyclonal antibody (CABT-BL2249)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptides derived from internal parts of Measles Virus Matrix protein.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Virus
<b>Purification</b>	Whole antiserum
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA, ICC
<b>Cellular Localization</b>	Envelope protein. Viral matrix protein. Virion
<b>Format</b>	Liquid
<b>Size</b>	200 µl
<b>Buffer</b>	Whole serum
<b>Preservative</b>	None
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

## GENE INFORMATION

Entrez Gene ID [1489803](#)

Protein Refseq [NP\\_056921](#)

---

UniProt ID [Q786M4](#)

---