



Anti-MAGEC1 (C-terminal) polyclonal antibody (CABT-BL2254)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide (conjugated with KLH) corresponding to C-terminus of human MAGEC1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Protein G purification
Conjugate	Unconjugated
Applications	WB, IHC, ELISA
Format	Liquid
Buffer	In PBS (0.09% sodium azide)
Preservative	0.09% Sodium Azide
Storage	Store at 4°C. For long term storage store at -20°C. Aliquot to avoid repeated freezing and thawing.

BACKGROUND

Introduction	This gene is a member of the melanoma antigen gene (MAGE) family. The proteins of this family are tumor-specific antigens that can be recognized by autologous cytolytic T lymphocytes. This protein contains a large number of unique short repetitive sequences in front of the MAGE-homologous sequence, and therefore is about 800 aa longer than the other
--------------	---

GENE INFORMATION

Entrez Gene ID	9947
----------------	----------------------

Protein Refseq	NP_005453
----------------	---------------------------

UniProt ID	O60732
------------	------------------------
