



Anti-MAGIX (aa 269-298) polyclonal antibody (CABT-BL2255)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, from within C terminal amino acids 269-298 of Human MAGIX (NP_079135.2).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Buffer	PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Entrez Gene ID	79917
Protein Refseq	NP_001093152
UniProt ID	A0A087WUY6