



Anti-MAL4P1.227 polyclonal antibody (CABT-BL2259)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fragment: IIHEHQEMEE TTEEGASIIIE RHSRTSSQNN DAGTSDYMYR ILSTTMEFPN EKETFGLYSP VSQTFDRALN FLRSVFISS SPVSQSKGQS KESQLAQIST(Plasmodium falciparum)
Isotype	IgG
Source/Host	Mouse
Species Reactivity	P. falciparum
Purification	Whole antiserum
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	50% Glycerol, Whole serum
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	812428
UniProt ID	Q8IFM2