



Anti-MAMDC4 (internal region) polyclonal antibody (CABT-BL2260)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	This antibody is specific for the Middle Region of the target protein.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA, IHC-P
Buffer	20 mM Potassium Phosphate, 150 mM Sodium Chloride, pH 7.0
Preservative	None
Storage	Store at -20°C. Avoid freeze-thaw cycles.

GENE INFORMATION

Entrez Gene ID	158056
Protein Refseq	NP_996803.2
UniProt ID	Q6UXC1