



Anti-MAML2 polyclonal antibody (CABT-BL2261)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH derived from within residues 1100 to the C-terminus of Human MAML2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ICC/IF
Cellular Localization	Nucleus speckle. Nuclear, in a punctate manner.
Format	Liquid
Buffer	1% BSA, PBS, pH 7.4
Preservative	0.02% Sodium Azide
Storage	Store at +4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

BACKGROUND

Introduction	Acts as a transcriptional coactivator for NOTCH proteins. Has been shown to amplify NOTCH-induced transcription of HES1. Potentiates activation by NOTCH3 and NOTCH4 more
---------------------	---

efficiently than MAML1 or MAML3.

GENE INFORMATION

Entrez Gene ID	84441
----------------	-----------------------

Protein Refseq	NP_115803
----------------	---------------------------

UniProt ID	Q8IZL2
------------	------------------------
