



# Anti-MAMSTR (aa 16-45) polyclonal antibody (CABT-BL2264)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                              |                                                                                                                     |
|------------------------------|---------------------------------------------------------------------------------------------------------------------|
| <b>Immunogen</b>             | Synthetic peptide conjugated to KLH, corresponding to a region within N terminal amino acids 16-45 of Human MAMSTR. |
| <b>Isotype</b>               | IgG                                                                                                                 |
| <b>Source/Host</b>           | Rabbit                                                                                                              |
| <b>Species Reactivity</b>    | Human                                                                                                               |
| <b>Purification</b>          | Immunogen affinity purified                                                                                         |
| <b>Conjugate</b>             | Unconjugated                                                                                                        |
| <b>Applications</b>          | WB                                                                                                                  |
| <b>Cellular Localization</b> | Nucleus.                                                                                                            |
| <b>Format</b>                | Liquid                                                                                                              |
| <b>Buffer</b>                | PBS                                                                                                                 |
| <b>Preservative</b>          | 0.09% Sodium Azide                                                                                                  |
| <b>Storage</b>               | Store at 4°C (up to 6 months). For long term storage store at -20°C                                                 |

## GENE INFORMATION

|                       |                           |
|-----------------------|---------------------------|
| <b>Entrez Gene ID</b> | <a href="#">284358</a>    |
| <b>Protein Refseq</b> | <a href="#">NP_872380</a> |

UniProt ID

[Q6ZN01](#)

---