



Anti-MAN2A2 (aa 113-162) polyclonal antibody (CABT-BL2267)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide from a region within N terminal amino acids 113-162 (PQDCQFALGG RGQKPELQML TVSEDLPF DN VEGGVWRQGF DISYSPNDWD) from Mouse MAN2A2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Golgi apparatus membrane.
Format	Liquid
Buffer	98% PBS, 2% Sucrose
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	140481
Protein Refseq	NP_766491
UniProt ID	Q197W7