



Anti-MAP6D1 polyclonal antibody (CABT-BL2291)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide from mouse SL21.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Affinity Purified
Conjugate	Unconjugated
Applications	WB
Preservative	None
Storage	6 months at 4°C from date of shipment

GENE INFORMATION

Entrez Gene ID	208158
Protein Refseq	NP_941001
UniProt ID	Q14BB9