



Anti-MAPT polyclonal antibody (CABT-BL2311)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	WB
Preservative	None
Storage	Store for 1 year at -20μ from date of shipment

BACKGROUND

Introduction	This gene encodes the microtubule-associated protein tau (MAPT) whose transcript undergoes complex, regulated alternative splicing, giving rise to several mRNA species. MAPT transcripts are differentially expressed in the nervous system, depending on stage of neuronal maturation and neuron type. MAPT gene mutations have been associated with several neurodegenerative disorders such as Alzheimers disease, Picks disease, frontotemporal dementia, cortico-basal degeneration and progressive supranuclear palsy. [provided by RefSeq, Jul 2008]
--------------	--

GENE INFORMATION

Entrez Gene ID	4137
Protein Refseq	NP_001116538

UniProt ID

[P10636](#)
