



# Anti-MARC2 (C-terminal) polyclonal antibody (CABT-BL2315)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	A KLH conjugated synthetic peptide selected from the C-terminal region of Human MOSC2.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB, ELISA
<b>Cellular Localization</b>	Mitochondrion outer membrane.
<b>Format</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Store at 4°C (up to 6 months). For long term storage store at -20°C

## BACKGROUND

<b>Introduction</b>	MOSC domain-containing protein 2, mitochondrial is a protein that in humans is encoded by the MOSC2 gene.
---------------------	---

# GENE INFORMATION

Entrez Gene ID [54996](#)

---

Protein Refseq [NP\\_060368](#)

---

UniProt ID [Q969Z3](#)

---