



Anti-MARC2 polyclonal antibody (CABT-BL2317)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within C terminal amino acids 303-333 of Human MOSC2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Cellular Localization	Mitochondrion outer membrane.
Format	Liquid
Buffer	99% PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

BACKGROUND

Introduction	MOSC domain-containing protein 2, mitochondrial is a protein that in humans is encoded by the MOSC2 gene.
---------------------	---

GENE INFORMATION

Entrez Gene ID [54996](#)

Protein Refseq [NP_060368](#)

UniProt ID [Q969Z3](#)
