



Anti-MB21D2 (aa 441-490) polyclonal antibody (CABT-BL2326)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 441-490 (QSDGGDPNQP DDRLAKKLQQ LVTENPGKSI SVFINPDDVT RPHFRIDDKF) of Human C3orf59.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Buffer	2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Entrez Gene ID	151963
Protein Refseq	NP_848591

UniProt ID

[Q8IYB1](#)
